

S-Monovette® GlucoEXACT


Information for Physicians



Reliable diagnosis of diabetes and gestational diabetes

- Conforms to DDG*/DGGG** guidelines on gestational diabetes
- Ultra-fast glycolysis inhibition with fluid preparation
- Optimal glucose stabilization up to 96 h
- Protects against misdiagnosis
- Citrate/Fluoride is recommended for diabetes screening

S-Monovette® GlucoEXACT



05.1074.001

04.1945.001

The S-Monovette® GlucoEXACT can be **reliably stabilized up to an impressive 96 h** at room temperature. The citrate fluoride solution in the S-Monovette® mixes directly with the blood during the blood collection process. The mixing process is usually completed by inverting the tube three times. Precise filling of the S-Monovette® GlucoEXACT guarantees reliable sample quality and accurate pre-analytics. In order to obtain an accurate result, the S-Monovette® GlucoEXACT must be filled precisely.

Plasma glucose level and multiplication factor:

The plasma glucose value (raw glucose value) determined must be **multiplied by a factor of 1.16** due to blood thinning caused by the fluid preparation. The multiplication factor is defined using the volume ratio of blood volume to preparation volume.

Centrifugation

Standard preparation	2000 x g	2500 x g	3000 x g	3500 x g	4000 x g
S-Monovette® GlucoEXACT, (Citrate/Fluoride)	9 min	8 min	7 min	6 min	5 min

For further information, go to:

www.sarstedt.com/service/zentrifugation

Ordering Information

Order number	Description	Length x diameter	Packaging
05.1074.001	S-Monovette® GlucoEXACT 3.1 ml	66 x 11 mm	50/500
04.1945.001	S-Monovette® GlucoEXACT 3.1 ml	75 x 13 mm	50/500

* Deutsche Diabetes Gesellschaft (German Diabetes Society)

** Deutsche Gesellschaft für Gynäkologie und Geburtshilfe (German Society for Gynecology and Obstetrics)

¹ Kleinwechter et al Gestationsdiabetes mellitus (GDM). Evidenzbasierte Leitlinie zur Diagnostik - Therapie und Nachsorge, DDG 08.2011

² S3-Leitlinie Gestationsdiabetes mellitus (GDM), Diagnostik, Therapie und Nachsorge, 2. Auflage, DDG, DGGG-AGG 2018

³ Nauck et al Definition, Klassifikation und Diagnostik des Diabetes mellitus Diabetologie 2017; 12 (Suppl 2):S94-S100 (DDG Praxisempfehlung)

⁴ Petersmann et al Definition, classification and diagnostics of diabetes mellitus JLabMed 2018 43(3): 73-79

⁵ Sacks et al Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus, Diabetes Care 34, e61-e99 06.2011

⁶ Yagmur et al Effective inhibition of glycolysis in venous whole blood and plasma samples JLabMed 2012 36:169-177

⁷ Will et al Whitepaper Sarstedt S-Monovette® GlucoEXACT A blood collection device for stabilizing glucose levels for 96 hours 2016

SARSTEDT, Inc.
1025 St. James Church Road
P.O. Box 468
Newton, NC 28658 - 0468
Toll free: +1 800 257 5101
Fax: +1 828 465 4003
customerservice.us@sarstedt.com
www.sarstedt.com