

Weißer Multiply®-µStrips

Optimiert für die qPCR

NEU!



Durch weiß eingefärbte Gefäße wird eine maximale Fluoreszenzreflektion bei der qPCR erreicht. Im Vergleich zu transparenten Produkten können mit weiß eingefärbten Verbrauchsmaterialien i.d.R. 10-fach höhere Fluoreszenzwerte gemessen werden.

Durch die optimierten Multiply®-µStrips ergeben sich beachtliche Vorteile:

- Stark erhöhte PCR-Sensitivität besonders bei niedrigen Template-Konzentrationen
- Leichter reproduzierbare Ergebnisse
- Vermeidung von „Crosstalk“ zwischen den einzelnen Wells
- Verbessertes Signal/Hintergrund Verhältnis
- Volumenreduktion möglich - somit können Reagenzien-Kosten deutlich (bis zu 40 %) reduziert werden

Unsere optimierten PRC-Ketten sind DNA-frei, DNase/RNase-frei sowie frei von PCR-Inhibitoren und erzielen durch dünnwandige Wells einen optimalen Temperaturtransfer.

Für einen perfekt abgestimmten effektiven Verschluss bieten wir hochtransparente Multiply® Lid Strips an.

Bestellinformation

Bestell Nummer	Bezeichnung	Verpackung
72.985.092	Multiply®-µStrip 0,2 ml, weiß	120 IK/480 UK
65.989.002	Multiply® Lid Strips, hochtransparent	120 IK/480 UK

SARSTEDT AG & Co.
Postfach 12 20 · D-51582 Nümbrecht
Telefon +49 22 93 305-0
Telefax +49 22 93 305-3450
☎ Service 0800 (Deutschland)
Telefon (0800) 0 83 305-0
info@sarstedt.com
www.sarstedt.com

 **SARSTEDT**

White Multiply[®]- μ Strips

Optimised for qPCR

NEW!



White tubes enable maximum fluorescence reflection in qPCR. Comparing to transparent products, 10-fold fluorescence levels can generally be measured with white disposables.

The optimised Multiply[®]- μ Strips provide considerable advantages and benefits:

- Convincingly improved PCR sensitivity especially at low template concentrations
- Enhanced reproducibility of results
- Prevention of cross talk between the wells
- Improved signal/background ratio
- Volume reduction is possible – for a significant decrease (up to 40%) in reagent costs

Our optimised PCR strips are DNA-free, DNase/RNase-free and PCR inhibitor-free, enabling optimal temperature transfer thanks to the thin-walled wells.

For effective closure, our range includes compatible, highly transparent Multiply[®] Lid Strips.

Ordering information

Order No.	Description	Packaging
72.985.092	Multiply [®] - μ Strip 0.2 ml, white	120 per box/480 per case
65.989.002	Multiply [®] Lid Strips, highly transparent	120 per box/480 per case

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 305-0
Fax +49 2293 305-3992
export@sarstedt.com
www.sarstedt.com

 **SARSTEDT**