

SMO

S-Monovette® Organiser

NEW!



Automated System for efficient, patient-specific preparation of blood collection

- Selects blood collection tubes corresponding to test requests
- Prints bar code labels
- Labels S-Monovettes
- Provides labels for additional tests
- Integrated bag packaging system ensures safe transfer to patient

Description and functions

- Stand-Alone System for the provision of patient-related, pre-labelled Sarstedt blood collection tubes (S-Monovettes)
- Label printing with barcode and text
- Tube labelling
- Provision of single labels for individual test requests
- Label size 40 x 26 mm
- Integral bag packaging system
- Offline operation with request for tubes via touch screen
- Online operation directly linked with LIS or HIS

Input:

- Tubes supplied in S-Monovette® Rack
- 50 tubes/rack
- 12 racks max.

Output:

- Transparent PE bags with 10 labelled tubes max., and separate labels
- Additional bags per order optionally available
- Tamper-evident closure
- Bags in easy reach for convenient handling
- Alternatively, bags may be unloaded on belt conveyor or desk/loading platform

Tube Types

- Outer diameter: 11 to 16 mm, length 65-100 mm (w/o cap)
- S-Monovettes

Dimensions

- 850 x 1424 x 1561 mm (WxDxH)
- Weight: approx. 350 kg

Mains Supply

- 115-230 V AC \pm 10%, 50/60 Hz, 16 A

General data

- Noise level: <70 dB(A)
- Adm. ambient temperature: + 15°C to + 40°C
- Max. relative humidity: 80%, non-condensing

Throughput

- Up to 750 tubes/h

Optional Accessories

- Conveyor belt
- Scanner for patient identification
- Zebra label printer

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 3050
Fax +49 2293 305- 2470
info@sarstedt.com
www.sarstedt.com