


Tube Stripper AST-2



AST-2 Basic model



- Motor stripper for blood bag tubes
- Mains or battery operation
- Versions with and without auto-stop feature
- Benchtop instrument with auto-stop feature and floor contact

Tube Stripper AST-2


AST-2 Basic model

Blood bag tube stripper
for easy, non-tiring operation




AST-2S with auto-stop feature

Stripper with automatic
motor stop close to blood bag



AST-2F Desk top instrument

Stripper with auto-stop feature,
fixture and floor contact operation



Technical Data

Model	AST-2	AST-2S	AST-2F
Order No.	87.7709.029	87.7709.025	87.7709.035
Right- or left-handed operation	✓	Please indicate, if required	-
Overload protection	✓	✓	✓
Forward and reverse operation	✓	✓	✓
Auto-stop feature	-	✓	✓
Floor contact	-	-	✓
Power supply	✓	✓	✓
Battery	✓	✓	✓
Weight	620 g	640 g	1,480 g
Motor	9V 20W	9V 20W	9V 20W
Dimensions including mains adapter (L x W x H in mm)	195 x 39 x 62	195 x 47 x 62	195 x 65 x 62
Transport case	✓	✓	-
Accessories			
Universal Mains adapter 100V - 240V AC 50-60 Hz	87.7709.033	87.7709.033	87.7709.033
and/or battery set NiMh 7.2V 700 mA	87.7709.031	87.7709.031	87.7709.031
Universal Charger for battery set 100V - 240V AC 50-60Hz	87.7709.032	87.7709.032	87.7709.032

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 3050
Fax +49 2293 305-122
info@sarstedt.com
www.sarstedt.com