

Transfusion devices



- For the transfusion of blood products
- Use of fine-pored sieve filters for effective elimination of micro aggregates
- High filter capacity and high flow rate due to a large filter surface

Transfusion devices



- 1 TRANSMED 10** – Integrated cascading filter 175 µm, 40 µm and 10 µm
- 2 TRANSMED 40** – Integrated cascading filter 175 µm and 40 µm
- 3 TRANSFLOW** – Integrated fine filter with a pore size of 40 µm
- 4 TRANSFLOW-N** – Integrated fine filter with a pore size of 40 µm
- 5 Transfusion device 200 µm** – Integrated fine filter with a pore size of 200 µm

Ordering information

Order no.	Description	Pck.
79.8200.000	1 TRANSMED 10	100
79.8200.100	2 TRANSMED 40	100
79.8200.440	3 TRANSFLOW	100
79.8200.490	4 TRANSFLOW-N	100
79.8200.500	5 Transfusion device 200 µm	300

Distributed by:

SARSTEDT AG & Co. KG
 P.O. Box 12 20
 D-51582 Nümbrecht
 Phone: +49 2293 305 0
 Fax: +49 2293 305 3992
 export@sarstedt.com
 www.sarstedt.com