

Transfusion devices



- For the transfusion of blood products by means of gravity
- Use of fine-pored sieve filters for effective elimination of micro aggregates
- Meets the requirements of EN ISO 1135-4





- 1 TRANSMED 10** – Integrated cascading filter 175 µm, 40 µm and 10 µm
- 2 TRANSMED 40** – Integrated cascading filter 175 µm and 40 µm
- 3 TRANSFLOW** – Integrated fine filter with a pore size of 40 µm
- 4 TRANSFLOW-N** – Integrated fine filter with a pore size of 40 µm

Ordering information

Order no.	Description	Pck.
79.8200.000	1 TRANSMED 10	100
79.8200.100	2 TRANSMED 40	100
79.8200.440	3 TRANSFLOW	100
79.8200.490	4 TRANSFLOW-N	100

Distributed by:

SARSTEDT AG & Co. KG
 P.O. Box 12 20
 D-51582 Nümbrecht
 Phone: +49 2293 305 0
 Fax: +49 2293 305 3992
 export@sarstedt.com
 www.sarstedt.com