

# REGANESTH® Spinal Needles Pencil-Point



Superior manufacturing quality and well-designed features

- Ergonomically designed needle hub
- CSF chamber with magnifying effect
- Atraumatic pencil-point tip

## The Needle Hub

- Ergonomic design for safe and convenient handling
- Magnifying effect enables the identification of even small fluid volumes



## The Pencil-Point Tip

- Atraumatic tip minimises injuries during puncture
- Reduces post-puncture headache
- Lateral rounded opening to prevent intra-arachnoidal malposition



## Ordering Information

Order No.	Description	Colour	Packaging
79.8081.520	REGANESTH® Spinal needle pencil-point 27 G x 38 mm	●	20
79.8081.500	REGANESTH® Spinal needle pencil-point 27 G x 90 mm	●	20
79.8081.515	REGANESTH® Spinal needle pencil-point 27 G x 90 mm with Introducer 22 G x 32 mm	●	20
79.8081.503	REGANESTH® Spinal needle pencil-point 27 G x 103 mm with Introducer 22 G x 32 mm	●	20
79.8081.510	REGANESTH® Spinal needle pencil-point 27 G x 110 mm	●	20
79.8081.530	REGANESTH® Spinal needle pencil-point 27 G x 120 mm	●	20
79.8081.540	REGANESTH® Spinal needle pencil-point 27 G x 120 mm with Introducer 20 G x 40 mm	●	20
79.8081.800	REGANESTH® Spinal needle pencil-point 25 G x 90 mm	●	20
79.8081.805	REGANESTH® Spinal needle pencil-point 25 G x 90 mm with Introducer 20 G x 30 mm	●	20
79.8081.803	REGANESTH® Spinal needle pencil-point 25 G x 103 mm with Introducer 20 G x 30 mm	●	20
79.8081.840	REGANESTH® Spinal needle pencil-point 25 G x 120 mm with Introducer 20 G x 40 mm	●	20
79.8081.811	REGANESTH® Spinal needle pencil-point 24 G x 90 mm with Introducer 20 G x 30 mm	●	20
79.8081.812	REGANESTH® Spinal needle pencil-point 24 G x 120 mm with Introducer 20 G x 40 mm	●	20
79.8081.210	REGANESTH® Spinal needle pencil-point 22 G x 90 mm	●	20

Now also available with the NRFit™ Non-Luer Connector in accordance with ISO 80369-6.



**NRFit™**  
Non-Luer Connector in accordance  
with ISO 80369-6

Distributed by:

SARSTEDT AG & Co. KG  
P.O. Box 12 20  
D-51582 Nümbrecht  
Phone: +49 2293 305 0  
Fax: +49 2293 305 3992  
export@sarstedt.com  
www.sarstedt.com