

# Storage Boxes

Convenient sample storage up to -90° C



For the convenient storage of Micro Tubes or Screw Cap Micro Tubes Sarstedt offers high quality storage boxes made from polypropylene. The robust boxes are available in six colours and ideally suited for long-time storage of 1.5 ml and 2.0 ml tubes with a maximum size of up to 40 x 12 mm.

## Advantages of the Storage Boxes at a glance

- Hinged lid or slip-on lid – convenient handling under any storage condition
- Alphanumeric coding at each aperture – facilitated sample identification
- Autoclavable – easy hygienic cleaning (121° C, 20 min.)
- Low-temperature-resistant – reliable storage up to -90° C
- Six colours – quick identification of the required storage box







### Slip-on lid

- Box size: 130 mm x 130 mm x 45 mm (L x W x H)
- Storage capacity 81 tubes
- Format 9 x 9

### Hinged lid

- Box size: 150 mm x 140 mm x 55 mm (L x W x H)
- Storage capacity 100 tubes
- Format 10 x 10

## Ordering Information

Box colour	Slip-on lid	Hinged lid	Packaging bag / case
Neutral 	93.876	93.877	5 / 20
Pink 	93.876.281	93.877.210	
Orange 	93.876.381	93.877.310	
Yellow 	93.876.481	93.877.410	
Green 	93.876.581	93.877.510	
Blue 	93.876.681	93.877.610	

SARSTEDT AG & Co. KG  
P.O. Box 12 20  
D-51582 Nümbrecht  
Phone +49 2293 305-0  
Fax +49 2293 305-3992  
export@sarstedt.com  
www.sarstedt.com