

Multi Sealer MS 500



- RF Hand-held sealer
- Safe sealing
- Power cord or battery operation



Multi Sealer MS 500

The new hand-held tube sealers, MS 500 M (mains) and MS 500 B (battery operated) are specially designed for use during blood donation or for the production of blood components.

The hand-held sealer features a specialised mechanical system that compensates for differences in user handling in order to ensure secure and consistent sealing.

The sealing head protector and electrode are also easy to detach for cleaning.

Product features

- Haemolysis-free sealing
- Precise seals
- Simple segment separation
- Indicator light and audible signal after sealing
- Easy electrode replacement
- Easy to clean
- Replaceable battery (MS 500 B)
- Durable hand sealer
- Compact case

Technical data

	MS 500 M (mains)	MS 500 B (battery operated)
Power supply	100 - 240 VAC	-
NiMH battery	-	24 V / 1700 mAh (replaceable)
Seals/charge	-	> 500
Full charge	-	abt. 3 hrs.
Charger	-	100 - 240 VAC, 41 VDC / 900 mA
Instrument dimensions (WxHxD)	18.0 x 5.5 x 23.0 cm	18.0 x 5.5 x 17.0 cm
Instrument weight	2.0 kg	2.0 kg
Hand sealer dimensions (Diameter x L)	3.1 x 18.5 cm	3.1 x 18.5 cm
Hand sealer weight	0.275 kg	0.275 kg
Case dimensions (WxHxD)	39.0 x 87.0 x 32.5 cm	39.0 x 87.0 x 32.5 cm
RF cable length	1.8 m	1.8 m

Ordering information

Order No.	Description
87.7708.850	MS 500 M (mains), includes hand-held sealer, RF cable, case
87.7708.860	MS 500 B (battery-operated), includes hand-held sealer, RF-cable, 1700 mAh battery, charger, strap, case

SARSTEDT AG & Co. KG
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 305 - 0
Fax +49 2293 305 - 3992
export@sarstedt.com
www.sarstedt.com

