



REGANESTH[®] unoplex

Ultrasound-guided peripheral nerve blocks




Suitable for ultrasound

In plane



Out of plane



REGANESTH[®] unoplex

- The REGANESTH[®] unoplex needle was developed for ultrasound-guided nerve blocks
- The needle tip enables outstanding visibility, even for steep puncture angles

Efficient

The *rhombus cut* ultrasound marker enables excellent echogenicity



Easy to use

Convenient handling thanks to the ergonomically designed REGANESTH® handle and bevel position marking



Accessories



Ordering information

Order no.	Description	Packaging unit
79.8096.050	REGANESTH® unoplex 22G x 50mm with detachable 60cm injection tube	10
79.8096.080	REGANESTH® unoplex 22G x 80mm with detachable 60cm injection tube	10
79.8096.100	REGANESTH® unoplex 20G x 100mm with detachable 60cm injection tube	10
79.8096.185	Unoplex catheter set 18G x 50mm/21G REGANESTH® unoplex 18G x 50mm Plexus catheter 21G x 400mm with stylet and hub insert, attachable 60cm injection tube, Tuohy-Borst adapter, flat filter 0.2 µm, REGAfix F flat filter fixation device	5
79.8096.150	REGANESTH® unoplex 18G x 50mm	20
79.8096.188	Unoplex catheter set 18G x 80mm/21G REGANESTH® unoplex 18G x 80mm Plexus catheter 21G x 400mm with stylet and hub insert, attachable 60cm injection tube, Tuohy-Borst adapter, flat filter 0.2 µm, REGAfix F flat filter fixation device	5
79.8096.180	REGANESTH® unoplex 18G x 80mm	20

Distributed by:

SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 305-0
Fax +49 2293 305-3992
export@sarstedt.com
www.sarstedt.com