

5 ml Micro Tubes

The intermediate volume option for labs

NEW!



The Micro Tube for medium volumes:

- Convenient handling of sample volumes between 2 and 5 ml
- Reliable SafeSeal locking closure
- PCR Performance quality for high purity requirements



SARSTEDT

5 ml Micro Tubes

Reliable handling of filling volumes between 2 and 5 ml

Routine applications in the lab like dilution series, cell breakdown or sample backup frequently require filling volumes above 2 ml. Due to their design or screw caps, larger volume tubes are sometimes not ideally suited for these purposes.

The new Sarstedt 5 ml Micro Tube combines the advantages of the SafeSeal Micro Tubes and 5 ml screw cap tubes in a single product. Easy to open with just one hand thanks to the ergonomic design, the snap-on lids with lip seal and retaining cams provide safe sealing characteristics even under thermal strain. The wide hinge enables easy and precise closure, and cap mix-ups associated with screw caps are eliminated.



Advantages of the 5 ml Micro Tube at a glance:

- Robust polypropylene – centrifugation to 21,000 x g, autoclavable
- Molded graduation – easy volume control
- Writing space and frosted lid – easy marking for quick identification
- Low height – minimizes contamination risk from pipette cone contact
- Compact size – space-saving backup storage of samples, media or cultures

Ordering Information

Order No.	Description	Packaging
72.701	5 ml Micro Tube	100 / bag 200 / case
72.701.400	5ml Micro Tube, PCR Performance Tested	100 / bag 200 / case

Accessories

Order No.	Description	Format	External dimensions (mm)	Packaging
93.850.001	Polycarbonate rack for 18 pcs. 5ml Micro Tubes	6 x 3	137 x 70 x 40 (L x W x H)	1 / case
93.1431.001	Polycarbonate rack for 48 pcs. 5ml Micro Tubes	12 x 4	257 x 90 x 40 (L x W x H)	1 / case

Sarstedt Inc.
1025 St. James Church Road
P.O. Box 468
NEWTON, NC 28658
Tel: +1 800 257 5101
Tel: +1 828 465 4000
Fax: +1 828 465 4003
customerservice@sarstedt.us
www.sarstedt.com