

# Sarstedt Cell Strainer

Quick and easy filtration of single-cell suspensions

**NEW!**



The Sarstedt cell strainer provides an innovative, fast and easy-to-handle solution for single-cell suspensions (primary cell cultures, flow cytometry).

## Features and Benefits

- Color-coded pore sizes – 40 µm (blue), 70 µm (white) and 100 µm (yellow)
- Circumferential rim with handle for easy and safe handling
- Four cams on the frame ensure continuous ventilation during filtration to prevent the ‘air lock’ effect and the cell strainer from overflowing
- Stackable design enables multi-stage filtration
- Compatible with standard 50 ml centrifuge tubes
- Adapter enables compatibility with a large number of smaller tubes (e.g. 15 ml centrifuge tubes, 12 x 75 mm FACS tubes)

## Ordering information

Order number	Description	Packaging
83.3945.040	Cell Strainer, 40 µm pore size, color code: blue, sterile	1/blister, 50/case
83.3945.070	Cell Strainer, 70 µm pore size, color code: white, sterile	1/blister, 50/case
83.3945.100	Cell Strainer, 100 µm pore size, color code: yellow, sterile	1/blister, 50/case
83.3945.999	Cell Strainer adapter for 15 ml, 5 ml, 17 x 100 mm and 12 x 75 tubes, sterile	1/blister, 25/case

Sarstedt, Inc.  
1025 St. James Church Road  
P.O. Box 468  
Newton, NC 28658 - 0468  
Tel: +1-828-465-4000  
Toll free: +1-800-257-5101  
Fax: +1-828-465-4003  
customerservice.us@sarstedt.com  
www.sarstedt.com

 **SARSTEDT**