

Safety-Multifly®-Kanüle ohne Multi-Adapter

Ideal für die Kurzzeitinfusion und Adaption an Luer-Verbindungen



Die Safety-Multifly®-Kanüle ohne den vormontierten Multi-Adapter kann direkt für die Kurzzeitinfusion sowie für die Konnektierung an Luer-Verbindungen eingesetzt werden.

Bestellinformation

Bestell Nummer	Bezeichnung	Ø / Länge	Ausführung	Farbcode	Schlauchlänge	Stk./ IK	Stk./ Karton
85.1637.201	Safety-Multifly®-Kanüle 20G lang ohne Adapter	0,9 mm / 19 mm	20G x 3/4", Nr. 1	gelb	240 mm	120	480
85.1637.203	Safety-Multifly®-Kanüle 20G kurz ohne Adapter	0,9 mm / 19 mm	20G x 3/4", Nr. 1	gelb	80 mm	120	480
85.1638.201	Safety-Multifly®-Kanüle 21G lang ohne Adapter	0,8 mm / 19 mm	21G x 3/4", Nr. 2	grün	240 mm	120	480
85.1638.203	Safety-Multifly®-Kanüle 21G kurz ohne Adapter	0,8 mm / 19 mm	21G x 3/4", Nr. 2	grün	80 mm	120	480
85.1640.201	Safety-Multifly®-Kanüle 23G lang ohne Adapter	0,6 mm / 19 mm	23G x 3/4", Nr. 16	blau	240 mm	120	480
85.1640.203	Safety-Multifly®-Kanüle 23G kurz ohne Adapter	0,6 mm / 19 mm	23G x 3/4", Nr. 16	blau	80 mm	120	480
85.1642.201	Safety-Multifly®-Kanüle 25G lang ohne Adapter	0,5 mm / 19 mm	25G x 3/4", Nr. 18	orange	240 mm	120	480

Safety-Multifly® Needle without Multi-Adapter

Ideal for short-term infusion and adaptation to Luer connection



The Safety-Multifly® Needle without the pre-assembled Multi-Adapter can be directly used for short-term infusion as well as for connection to Luer connections.

Ordering information

Order No.	Description	Ø / Length	Version	Colour code	Tube length	Pcs./ inner case	Pcs./ case
85.1637.201	Safety-Multifly®-Needle 20G long without adapter	0.9 mm / 19 mm	20G x ¾", No. 1	yellow	240 mm	120	480
85.1637.203	Safety-Multifly®-Needle 20G short without adapter	0.9 mm / 19 mm	20G x ¾", No. 1	yellow	80 mm	120	480
85.1638.201	Safety-Multifly®-Needle 21G long without adapter	0.8 mm / 19 mm	21G x ¾", No. 2	green	240 mm	120	480
85.1638.203	Safety-Multifly®-Needle 21G short without adapter	0.8 mm / 19 mm	21G x ¾", No. 2	green	80 mm	120	480
85.1640.201	Safety-Multifly®-Needle 23G long without adapter	0.6 mm / 19 mm	23G x ¾", No. 16	blue	240 mm	120	480
85.1640.203	Safety-Multifly®-Needle 23G short without adapter	0.6 mm / 19 mm	23G x ¾", No. 16	blue	80 mm	120	480
85.1642.201	Safety-Multifly®-Needle 25G long without adapter	0.5 mm / 19 mm	25G x ¾", No. 18	orange	240 mm	120	480

Technical modifications reserved

SARSTEDT AG & Co. KG
P.O. Box 12 20
D-51582 Nümbrecht
Phone: +49 2293 305 0
Fax: +49 2293 305 3992
export@sarstedt.com
www.sarstedt.com

 **SARSTEDT**