

# 兼容性表 - 移液吸头

## 带滤芯的移液吸头

容量	产品编号	兼容性
2.5 µl	70.3010.265, 70.3010.285, 70.3010.365, 70.3010.385	● 2.5, ● 2.5
10 µl	70.3010.255, 70.3010.275, 70.3010.355, 70.3010.375	● 10
20 µl	70.3020.255, 70.3020.275, 70.3020.355, 70.3020.375	● 20
20 µl	70.3030.265, 70.3030.285, 70.3030.365, 70.3030.385	● 20, ● 20
100 µl	70.3030.255, 70.3030.275, 70.3030.355, 70.3030.375	● 100
200 µl	70.1189.215	● 200
200 µl	70.3031.255, 70.3031.275, 70.3031.355, 70.3031.375	● 200
300 µl	70.3040.255, 70.3040.275, 70.3040.355, 70.3040.375	● 300
1,000 µl	70.3050.255, 70.3050.275, 70.3050.355, 70.3050.375	● 1,000
1,000 µl	70.3060.255, 70.3060.275, 70.3060.355, 70.3060.375	● 1,000
1,250 µl	70.1186.210, 70.1186.410	● 1,250

## 无滤芯的移液吸头

容量	产品编号	兼容性
10 µl	70.3010, 70.3010.200, 70.3010.210, 70.3010.205, 70.3010.100, 70.3010.110	● 10
20 µl	70.3020, 70.3020.200, 70.3020.210, 70.3020.205, 70.3020.100, 70.3020.110	● 20
20 µl	70.3021, 70.3021.200, 70.3021.210, 70.3021.205, 70.3021.100, 70.3021.110	● 20
200 µl	70.3030, 70.3030.020, 70.3030.210, 70.3030.205, 70.3030.100, 70.3030.110	● 200
200 µl	70.1189.105	● 200
250 µl	70.3031, 70.3031.020, 70.3031.210, 70.3031.205, 70.3031.100, 70.3031.110	● 250
300 µl	70.3040, 70.3040.100, 70.3040.210, 70.3040.110, 70.3040.205, 70.3040.305	● 300
1,000 µl	70.3050, 70.3050.020, 70.3050.210, 70.3050.030, 70.3050.100, 70.3050.110, 70.3050.205, 70.3050.305	● 1,000
1,250 µl	70.3060, 70.3060.100, 70.3060.110, 70.3060.210, 70.3060.205, 70.3060.305	● 1,250
1,250 µl	70.1186, 70.1186.100, 70.1186.200	● 1,250

- = 兼容
- 10 = 移液容积有限
- ✗ = 不推荐
- = 未进行验证
- MC = 多通道移液器

移液器制造商会根据技术条件的变化来更改移液器顶锥。因此，我们在兼容性表中所提出的建议仅供参考。如果您想检查我们的吸头与您所使用的移液器的匹配度以及卸载性能，我们可免费为您提供吸头样品。

## Eppendorf® 移液器

容量	Eppendorf Reference®	Eppendorf Research®	Eppendorf Research® plus	Eppendorf Research® pro	Xpert®
0.1-2.5 µl	● 2.5	● 2.5	● 2.5	● 2.5	
0.5-10 µl	● 10	✗	● 10	● 10	
2-20 µl	● 20	● 20	● 20	● 20	
10-100 µl		✗ 20	● 20	● 20	
20-200 µl		● 20	● 20	● 20	
30-300 µl MC		● 300	● 300	● 300	
100-1,000 µl		● 100	● 100	● 100	
500-5,000 µl		● 500	● 500	● 500	
1-10 ml					
0.1-2.5 µl MC					
0.5-10 µl MC					
2-20 µl MC					
10-100 µl MC					
20-200 µl MC					
30-300 µl MC					
100-1,000 µl MC					
500-5,000 µl MC					
0.5-10 µl MC					

## Thermo® 微量移液器

容量	Finnpipette® Digital / F1 / F2	Finnpipette® Focus	Finnpipette® Novus	Bio Control	Pipetman® / Ultra / Concept
0.2-2 µl	● 2.5	● 2.5			● 2.5
F2 0.5-5 µl	● 5				✗
F1 1-10 µl	● 10				✗
0.5-10 µl Univ.	● 10				
0.5-10 µl MC	● 10				
2-20 µl Univ.	● 20				
5-50 µl MC	● 50				
10-100 µl Univ.	● 100				
20-200 µl Univ.	● 200				
50-300 µl MC	● 300				
100-1,000 µl	● 1,000				
200-1,000 µl	● 2,000				
0.3-3 µl		● 2.5			
0.5-5 µl		● 5			
1-10 µl		● 10			
30-300 µl		● 300			
100-1,000 µl		● 1,000			
30-300 µl MC		● 300			
30-300 µl MC		● 300			
100-1,000 µl MC		● 1,000			
100-1,200 µl MC		● 1,200			
50-300 µl MC		● 500			
P20/U2 0.2-2 µl					● 2.5
P10/U10 0.1-10 µl					● 10
C8/C12 1-10 µl					● 10
P20/U20 2-20 µl					● 20
P100/U100 20-100 µl					● 100
P200/U200 20-200 µl					● 200
U8 20-300 µl					● 200
P1,000/U3,000 200-1,000 µl					● 1,000
P5,000/U5,000 1,000-5,000 µl					● 5,000
P10,010 1-10 ml					● 10

