


ピペットチップ適合表


フィルターチップ



70.1130.212 70.1130.217
70.1130.210 70.1130.215
70.1114.210 70.1114.215
70.1116.210 70.1116.215
70.760.213 70.760.219
70.760.212 70.760.217
70.1189.215
70.760.211 70.760.216
70.3040.255 70.3040.275
70.3050.255 70.3050.275
70.3060.255 70.3060.275
70.1186.210 70.1186.410

Eppendorf® Pipette				Thermo® Finnpipette							Gilson® Pipette						
Eppendorf Reference®		Eppendorf Research®		Eppendorf Research® plus		Eppendorf Research® pro		Finnpipette® Digital / F1 / F2		Finnpipette® Focus		Finnpipette® Novus		Pipetman® / Ultra / Concept			
● 2.5	● 2.5	● 2.5		● 2.5		● 2.5	● 2.5	●	● 2.5	● 2.5				● 2.5	×		
● 10		×	●	×	●	×	●	●	●	●				●	×		
● 20	● 20	● 20		● 20	● 20	● 20	● 20	●	●	●				●	×		
● 20	● 20	● 20		● 20	● 20	● 20	● 20	●	●	●				●	×		
● 100	● 100	● 100		● 100	● 100	● 100	● 100	●	●	●				●	×	●	
● 200	● 200	● 200		● 200	● 200	● 200	● 200	●	●	●				●	×	●	
● 300	● 300	● 300		● 300	● 300	● 300	● 300	●	●	●				●	×	●	
● 1,000	● 1,000	● 1,000		● 1,000	● 1,000	● 1,000	● 1,000	●	●	●				●	×	●	
● 1,250	● 1,250	● 1,250		● 1,250	● 1,250	● 1,250	● 1,250	●	●	●				●	×	●	

ピペットチップ



70.3010 70.1130.100 70.1130.105
70.3020 70.1114.100 70.1114.200
70.3021 70.1116.100
70.3030 70.760.102 70.760.107
70.1189.105
70.3031 70.760.201
70.3040 70.3040.100 70.3040.110 70.3040.200
70.3050 70.3050.020 70.3050.030 70.3050.100 70.3050.110
70.3060 70.3060.100 70.3060.110 70.3060.200
70.1186 70.1186.100 70.1186.200

Eppendorf Reference®				Eppendorf Research®		Eppendorf Research® plus		Eppendorf Research® pro		Finnpipette® Digital / F1 / F2		Finnpipette® Focus		Finnpipette® Novus		Pipetman® / Ultra / Concept			
● 10	● 10	● 10		● 10		● 10	● 10	●	●	●				●	×				
● 20	● 20	● 20		● 20	×	● 20	● 20	●	●	●				●	×				
● 20	● 20	● 20		● 20	● 20	● 20	● 20	●	●	●				●	×				
● 200	● 200	● 200		● 200	● 200	● 200	● 200	●	●	●				●	×	●			
● 200	● 200	● 200		● 200	● 200	● 200	● 200	●	●	●				●	×	●			
● 250	● 250	● 250		● 250	● 250	● 250	● 250	●	●	●				●	×	●			
● 300	● 300	● 300		● 300	● 300	● 300	● 300	●	●	●				●	×	●			
● 1,000	● 1,000	● 1,000		● 1,000	● 1,000	● 1,000	● 1,000	●	●	●				●	×	●			
● 1,250	● 1,250	● 1,250		● 1,250	● 1,250	● 1,250	● 1,250	●	●	●				●	×	●			

- = 適合
- 10 = 記載容量まで適合
- ×
- MC = マルチチャンネルピペット
- = データ無

適合表に記載されている情報は推奨されている一例であり、各メーカーの仕様変更等により適合・不適合が変更されている場合があります。適合確認が必要な場合には実際にお試しください。

