

# Multi-Safe vario 2000

Large opening, 2 litre



## The innovative solu


With innovative day and end closure through one-handed operation!



Clear fill level marking thanks to new label design



Compatible with standard trays




77.3893.020

Disposal opening can be... by breaking ou

Luer systems  
Needles

Easy handling during lid assembly due to snap-on hooks



# Option for safe disposal

Integrated QR code for direct online access to operating instructions




Can also be ordered with Biohazard labelling!



77.3893.120




individually adjusted  
cut sections



Scalpel  
Luer systems

Integrated stamp in the lid allows easy and safe uncoupling of the insulin pen when closing the lid




## Ordering information

Order no	Description	Volume	Material	Pieces per case	Units per pallet	Case dimensions	Pallet dimensions
77.3893.020	Multi-Safe vario 2000	2.0 l	PP	50	800	59.5 x 39.4 x 38.8 cm	120 x 80 x 195 cm
77.3893.120	Multi-Safe vario 2000 with biohazard marking	2.0 l	PP	50	800	59.5 x 39.4 x 38.8 cm	120 x 80 x 195 cm

- Individually extendible and washable brackets

Wall bracket



Round tube holder



Fastening plate



## Ordering information for accessories

Order no	Description	Size	Material	Pieces per case	Case dimensions
95.963.020	Wall holder vario 2000	23.5 x 10.8 x 5.8 cm	Stainless steel	5	26 x 19 x 9.5 cm
95.963.006	Round tube holder for wall holder	13.0 x 4.0 x 9.0 cm	POM	10	28 x 21 x 13.5 cm
95.963.005	Fastening plate for wall holder	13.0 x 7.9 x 1.3 cm	Stainless steel	20	26 x 18.3 x 8.5 cm



Further information on our extensive range of disposal boxes and the complete brochure 412 can be found at [www.sarstedt.com](http://www.sarstedt.com).

SARSTEDT AG & Co. KG  
P.O. Box 12 20  
D-51582 Nümbrecht  
Phone: +49 2293 305 0  
Fax: +49 2293 305 3992  
[export@sarstedt.com](mailto:export@sarstedt.com)  
[www.sarstedt.com](http://www.sarstedt.com)