

Blood Collection in Paediatrics



Blood collection systems for any age

- S-Monovette® & Safety-Multifly®-Needle
- Microvette® & Safety-Lancet
- Micro Tubes & Micro-Needle



Blood collection systems for any age

Blood collection in paediatrics calls for skilled phlebotomists and high-standard collection systems. The Sarstedt range includes a wide variety of different collection systems for venous and capillary blood collection from neonates and young children.

S-Monovette[®], Safety-Needle & Safety-Multifly[®]-Needle

As the sensitivity of modern analyzer systems increases, the adequate sample volume can be reduced to a minimum also for extensive routine tests.

With its special design of reduced dimensions and low nominal volume the S-Monovette[®] Paediatrics is ideally suited to accommodate these sophisticated requirements.

The aspiration principle permits gentle blood collection so that the S-Monovette[®] Paediatrics fully complies with requirements imposed by critical vein conditions.



Microvette[®] & Safety-Lancet

Depending on the application, the Microvette[®] is available with a cylindrical or conical insert and a range of volumes from 100 to 500 µl. Blood can be collected either using capillary tubes, or by venipuncture using the tube rim.


The Safety-Lancet ensures optimum safety for both the patient and user and is available in a range of sizes either with a needle or blade.



Micro Tubes & Micro-Needle

The Micro-Needle has been designed to accommodate the special requirements in neonatology. It features a ridged grip for optimal and secure needle control. The short needle reduces blood waste to a minimum.

Compatible sample tubes in a range of preparations with screw cap or push cap are available.



SARSTEDT AG & Co.
P.O. Box 12 20
D-51582 Nümbrecht
Phone +49 2293 3050
Fax +49 2293 305-122
info@sarstedt.com
www.sarstedt.com